



US005698155A

# United States Patent [19]

Grosswald et al.

[11] Patent Number: 5,698,155

[45] Date of Patent: Dec. 16, 1997

- [54] **METHOD FOR THE MANUFACTURE OF PHARMACEUTICAL CELLULOSE CAPSULES**
- [75] Inventors: **Ralph R. Grosswald**, Fairfield, Iowa; **Jeffory B. Anderson**, Springville; **Clair S. Andrew**, Provo, both of Utah
- [73] Assignee: **GS Technologies, Inc.**, Fairfield, Iowa

- 97100247 5/1997 European Pat. Off. .
- 8027378 9/1981 Japan .
- 87120178 11/1988 Japan .
- 87120179 11/1988 Japan .
- 89145626 12/1990 Japan .
- 84222676 4/1996 Japan .
- WO91/19487 12/1991 WIPO .

- [21] Appl. No.: 377,669
- [22] Filed: Jan. 24, 1995

### Related U.S. Application Data

- [63] Continuation of Ser. No. 893,091, May 29, 1992, abandoned, which is a continuation-in-part of Ser. No. 708,023, May 31, 1991, abandoned.
- [51] Int. CL<sup>5</sup> ..... B29C 33/06; B29C 41/14; B29C 41/46
- [52] U.S. CL ..... 264/402; 264/405; 264/297.8; 264/305; 264/310; 264/327; 264/DIG. 37; 425/174.6; 425/272; 425/804
- [58] Field of Search ..... 425/269, 270, 425/272, 274, 275, 804, 174, 174.4, 174.6, 547, 577; 264/215, 299, 303, 304, 310-311, 327, DIG. 37, 402, 405

### [56] References Cited

#### U.S. PATENT DOCUMENTS

1,125,619	1/1915	Winchester	.....	425/96
1,787,777	1/1931	Colton	.....	425/804
1,978,829	10/1934	Wilkie	.....	264/303
2,254,262	9/1941	Bratring	.....	264/303
2,299,039	10/1942	Scherer	.....	
2,526,683	10/1950	Murphy	.....	264/304

(List continued on next page.)

#### FOREIGN PATENT DOCUMENTS

0056825	8/1982	European Pat. Off. .
0401832 A	12/1990	European Pat. Off. .

### OTHER PUBLICATIONS

- Parke, R.F., "Hydroxypropylmethylcellulose Capsules", Oct. 20, 1966.
- Ridgway, K., Ed. (1987) "Hard Capsules Development and Technology," The Pharmaceutical Press, London, pp. 56-58.

Primary Examiner—Angela Ortiz  
 Attorney, Agent, or Firm—John L. Lee

### [57] ABSTRACT

A method and apparatus for manufacturing pharmaceutical capsules use an aqueous solution of a thermogelling cellulose ether composition and use capsule body pins and capsule cap pins as molds. The method involves heating the pins, dipping the pins into the solution to cause the solution to gelatinize on the surface of the pins, removing the pins and drying the gelatinized solution on the surface of the pins to form capsule bodies and capsule caps. Pins are heated pre-dip and post-dip to facilitate gelating. Counterflow air is applied to provide drying from the inside. Capsule parts are removed by gripping. Capsule parts may have a thick wall and a stiffening ring.

26 Claims, 15 Drawing Sheets

